

CPAPs (CONTINUOUS POSITIVE AIRWAY PRESSURE)

Specific requirements apply to the carriage of battery operated CPAP equipment. When travelling on Qantas Airways, QantasLink, The Jetstar Group, National Jet Systems or Network Aviation the following applies:

- lithium batteries, battery packs and spare batteries **not exceeding 100Wh** – do not require a dangerous goods operator (airline) approval;
- lithium batteries, battery packs and spare batteries **exceeding 100Wh but not exceeding 160Wh** – require a dangerous goods operator (airline) approval and spare batteries are limited to two per passenger (contact dg@qantas.com.au for more information)
- all spare batteries and battery packs/powerbanks must be carried in carry-on baggage only

Prior to travel, passengers must contact Qantas Specific Needs (specialhandling@qantas.com.au), to obtain a medical clearance and if applicable, request the use of aircraft power.

BATTERY OPERATED		
WATT HOUR (Wh) RATING = VOLTAGE (V) x AMP HOURS (Ah)		
Brand	Model	Watt Hours
ResMed	AutoSet Spirit	95.04Wh
	S7 Lightweight	59.94Wh
	S9 AutoSet	97.5Wh
	S9 Elite	70Wh
	Stellar 100,150	23Wh
ResMed (Sullivan)	APD 2	95Wh
	APD 2S	95Wh
	III	95Wh
	PD 2E	95Wh
Respironics	BIPAP Pro 2	36Wh
	BIPAP Pro Synchrony	66Wh
	Remstar M Series	36Wh
	Remstar Plus	97.2Wh
Invacare	XPO2	77Wh
Medical Industries America	C1001	36Wh
Puritan Bennett	425	40Wh
	420E/G/S	25Wh
	420G	25Wh
	420S	25Wh
Schiller	Argus PRO LifeCare	44.72Wh

MAINS POWERED	
Brand	Model
Resmed	Autoset T
	S6 – Elite, Lightweight , Plus
	S7 - Elite
	S8 – Compact , Lightweight, Escape , Escape II
	Sullivan Naval- I , II (all) , III (all) , V (all)
	Remstar – Auto , Auto A-Flex, LX, Pro-flex, Choice, Choice LX, Plus C-Flex, Solo LX, Solo Plus LX
Transcend	Auto
	II
Sunrise Medical	9000 / 9001 series
Laerdal Heartstart	FR2+
	FR3
	FR3
Healthdyne	7300 , 7300H

BATTERY PACKS		
Brand	Model	Watt Hours
Transcend	P4	38Wh
	P8	75Wh
Roscoe	Volt	32.56Wh
Resmed	Power Station II	97Wh

This is not a comprehensive list.

POCs (PORTABLE OXYGEN CONCENTRATORS)

Specific requirements apply to the carriage of battery operated CPAP equipment. When travelling on Qantas Airways, QantasLink, The Jetstar Group, National Jet Systems or Network Aviation the following applies:

- lithium batteries, battery packs and spare batteries **not exceeding 100Wh** – do not require a dangerous goods operator (airline) approval;
- lithium batteries, battery packs and spare batteries **exceeding 100Wh but not exceeding 160Wh** – require a dangerous goods operator (airline) approval and spare batteries are limited to two per passenger (contact dg@qantas.com.au for more information)
- all spare batteries and battery packs must be carried in carry-on baggage only

Prior to travel, passengers must contact Qantas Specific Needs (specialhandling@qantas.com.au), to obtain a medical clearance and if applicable, request the use of aircraft power.

BATTERY OPERATED		
WATT HOUR (Wh) RATING = VOLTAGE (V) x AMP HOURS (Ah)		
Brand	Model	Watt Hours
AirSep Corporation	Freestyle	97.68Wh
	Lifestyle	36Wh
DeVibiss Healthcare/Sunrise Medical	iGO (306DS)	98Wh
InoGen	One G2	98Wh
	One G3	84Wh
iNOVA	Life Choice	76.96Wh
	Life Choice Activox Sport	76.96Wh
Oxlife	Independence	95Wh
Respironics	EverGo (900)	98Wh
	Eclipse (1000)	97.5Wh
	Eclipse 2 (1000A)	97.5Wh
	Eclipse 3 (1000B)	97.5Wh
Philips Respironics	SimplyGo	95Wh
	PD 2E	95Wh
SeQual	Eclipse	97.5Wh x 2
	Eclipse 2	97.5Wh x 2
	Eclipse 3	97.5Wh x 2
	Eclipse 5	97.5Wh x 2
	eQuinox	95Wh x 2
International Biophysics Corp	Life Choice	72Wh

This is not a comprehensive list.

MAINS POWERED	
Brand	Model
Invacare	XP02 (XPO130)

BATTERY PACKS		
Brand	Model	Watt Hours
Invacare	XP0110	77Wh